
Oxygenator with integrated arterial filter / Oxygénateur avec filtre artériel intégré / Oxygenator mit integriertem arteriellem Filter / Oxigenador con filtro arterial integrado / Oxigenador com filtro arterial integrado / Ossigenatore con filtro arterioso integrato / Oxygenator met geïntegreerde arteriële filter / Oxygenator med inbyggd artärfilter / Oxygenator med integreret arteriefilter / Oksygenator med integrert arteriefilter / Hapetin, jossa on sisäänrakennettu valtimosuodatin / Оксигенатор со встроенным артериальным фильтром / Oksygenator ze zintegrowanym filtrem tętniczym / Entegre arterielyl filtreli oksijenaattori / Oksigenator sa integrisanim arteriiskim filterom / Оксигенатор з інтегруванням артеріальним фільтром

CAPIOX[®] FX 25

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

Read all warnings, cautions and instructions carefully prior to use.

DESCRIPTION OF PRODUCTS

The CAPIOX FX25 is a membrane oxygenator with microporous polypropylene hollow fibers, to be used as an extracorporeal gas exchange device in which the blood flows on the outside of the fibers and with the ventilating gas flowing through the fibers.

The CAPIOX FX25 consists of a module for gas exchange with an integrated heat exchanger and Arterial Filter. The CAPIOX FX25 is available with (Code No.: CX*FX25RW, CX*FX25RE) or without (Code No.: CX*FX25W, CX*FX25E) a Hardshell Reservoir with an integrated cardiomy filter to provide a single-unit construction for ease of use. (CX*FX25RW, CX*FX25W and CX*FX25RE, CX*FX25E differ in the positioning of the arterial blood port of the oxygenator.)

The blood contacting surfaces are coated with **Xcoating**.

Xcoating is bio-compatible material that is applied on the blood-contacting surfaces of the device to reduce the adhesion of platelets to the device.

INDICATION FOR USE

The CAPIOX FX25 is intended to be used during open heart surgical procedures requiring cardiopulmonary bypass for periods up to 6 hours. The CAPIOX FX25 is for use with patients when required blood flow rate will not exceed 7L/min.

The CAPIOX FX25 Hardshell Reservoir is also intended for use in vacuum assisted venous drainage procedure, in post-operative chest drainage and autotransfusion procedures to aseptically return the blood to the patient for blood volume replacement.

The Integral Arterial Filter is intended to filtrate non-biologic particles and emboli and to facilitate gaseous emboli removal from the blood flowing through a cardiopulmonary bypass circuit.





DIRECTION OF BLOOD OUTLET PORT

The oxygenator is available with TWO OUTLET PORT CONFIGURATIONS - LEFT side (W) or RIGHT side (E).

Refer to the diagram below to determine which configuration is needed for your circuit.

TYPE	LEFT	RIGHT
Oxygenator with Hardshell Reservoir	code : CX*FX25RW	code : CX*FX25RE
Oxygenator without Hardshell Reservoir	code : CX*FX25W	code : CX*FX25E
Place waterports facing South as indicated in the illustration. The blood outlet port orientation is (W)est - LEFT side or (E)ast - RIGHT side. Note: this illustration is a top view of the oxygenator.	Blood inlet port Blood outlet port (Left) Gas outlet port Water ports	Blood inlet port Gas outlet port Blood outlet port (Right) Water ports

STRUCTURE

[REF] : CX*FX25RW


[REF] : CX*FX25RE

[REF] : CX*FX25W

[REF] : CX*FX25E

CX*FX25RW, CX*FX25W and CX*FX25RE, CX*FX25E differ in the positioning of the arterial blood outlet port of the oxygenator. All figures hereafter only refer to FX25RW.

STRUCTURE

WARNINGS

Describe serious adverse reactions and potential safety hazards, limitations in use imposed by them, and steps that should be taken if they occur.

WARNINGS

- CAPIOX FX25 is designed to operate at blood flow rates within the range of 0.5 to 7.0 L/min. Do not use any blood flow rate outside this range.
- Do not use solvents such as alcohol, ether, acetone, etc. These solvents may cause damage if used in or on the device.
- In order to prevent gaseous emboli from entering the blood phase, follow the instructions below.
 - DO NOT OBSTRUCT GAS OUTLET PORT. Avoid build up of excess pressure in the gas phase to prevent gaseous emboli entering the blood phase.
 - Pressure in the blood phase should always be higher than that in the gas phase to prevent gaseous emboli entering the blood phase.
 - The gas flow rate should not exceed 20 L/min. Excessive gas flow rate will bring about pressure increase in the gas phase, allowing gaseous emboli to enter the blood phase.
 - During recirculation, do not use pulsatile flow and do not stop the blood pump suddenly as these actions may cause gaseous emboli to enter the blood phase from the gas phase due to inertia force.
 - When CAPIOX FX25 Oxygenator module is used separately from the Hardshell Reservoir, set the module so that the upper end of the fibers is lower than the blood level in the venous reservoir. This prevents gaseous emboli from entering the blood phase from the gas phase.
 - To prevent gaseous emboli from entering the blood phase, make sure that the arterial pump flow rate always exceeds the flow rate of the cardioplegia line. The blood flow rate of the cardioplegia line should not exceed 1 L/min.
 - Do not exceed a 15°C (27°F) temperature difference between blood and water in heat exchanger to prevent gas dissolved in blood from forming bubbles.
- Total flow rate of the arterial line and any separate arterial lines must not exceed the flow rate at the oxygenator inlet port.
- Pressure at the blood inlet of the oxygenating module should not exceed 133 kPa (1,000 mmHg). Pressures greater than 133 kPa (1,000 mmHg) may cause leaks or damage to the device.
- Water pressure at the heat exchanger inlet should not exceed 28 PSI (2 kgf/cm²) (196 kPa). Pressures greater than 28 PSI (2 kgf/cm²) (196 kPa) may cause leaks or damage to the device.
- Adequate heparinization of the blood is required to prevent it from clotting in the system.

- When using the centrifugal pump on the arterial line, clamp the arterial line distal to the oxygenator (the patient's side) before stopping the pump. Improper clamping may cause back-flow of blood or migration of gaseous emboli into the blood side.
- Use caution when removing air during priming and perfusion. Excessive shock to the device – especially with hard objects – can cause damage to the device.
- Do not reduce heparin during circulation. Otherwise, blood clotting might occur.

Use CAPIOX FX25 Oxygenator with Hardshell Reservoir taking note of the additional warnings below.

- The blood flow into the cardiotomy filter should not exceed the rate of 5 L/min. Excessive blood flow rate may increase the pressure in the cardiotomy filter, resulting in back-flow into any solution or blood administration line connected to the Hardshell Reservoir.
- Yellow vent port cap need not be removed, since this port assures adequate venting with the cap attached. Do not close the vent port, as this may cause positive pressure in the Hardshell reservoir, resulting in back-flow into solution or blood administration lines connected to the Hardshell Reservoir.
- Before stopping the pump, release the pressure inside of the reservoir to the atmospheric pressure. Otherwise, venous blood might be drained from the patient.
- Negative pressure below -20 kPa (-150 mmHg) should not be applied to this reservoir as this may damage it.
- When the level of blood in the reservoir is below the minimum operating level, blood flow at a high flow rate into the cardiotomy filter may generate gaseous micro emboli (GME) which could migrate to the patient.
- Minimum operating volume in the reservoir is 200 mL. Set appropriate blood storage level, relative to venous flow rate, to prevent gaseous emboli passing to patient. (Refer to Fig. 11, page 10 "SAFETY MARGIN")
- Do not use this product for a period in excess of six hours. Excessive use for over six hours may lead to plasma leak and thrombi formation, which may compromise the gas exchange performance.
- Do not change the position of the positive pressure relief valve on the hardshell reservoir. Placing the positive pressure relief valve in any other location on the reservoir could expose the positive pressure relief valve to fluid that could impair the performance of the valve and could allow positive pressure build up to occur in the reservoir. If this occurs during VAVD, it may result in reduced venous drainage and retrograde flow of air into the venous line.

PRECAUTIONS

Include information regarding any special care to be exercised by the practitioner for the safe and effective use of the device.

PRECAUTIONS

- Read all the instructions of the drugs and the medical devices to be used in combination with this device prior to use.
- This device should only be used by properly trained and qualified personnel.
- This device is sterile and non-pyrogenic in the unopened, undamaged package. Inspect the device and package carefully. Do not use if the package and/or device is damaged, or if caps are not in place.
- Dispose of safely after single use to avoid risk of infection.
- This product has been sterilized by ethylene oxide gas.
- For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- In case of the need to replace the oxygenator, make sure that there is a spare oxygenator available.
- Use aseptic technique in all procedures.
- Perform priming using crystalloid solution which contains no blood, plasma and/or blood derived products.
- This device should be used with an appropriate holder supplied by TERUMO.
- Place CAPIOX FX25 Oxygenator module downstream from the blood pump.
- Thermistor probes can be used with a temperature monitor from Measurement Specialties Inc. by use of the cables supplied.
- Temperature of water delivered to heat exchanger must not exceed 42°C (108°F), otherwise blood damage may result.
- A minimum of 0.5 L/min oxygen gas flow is needed when blood is circulated. Less than 0.5 L/min oxygen gas flow may result in inadequate gas exchange.
- Stop the gas flow when circulation is suspended. During recirculation, check blood gas pressure. Excessive gas flow may cause low PaCO₂, alkalosis, or blood damage.
- Upon patient rewarming, adjust O₂ concentration, gas flow rate and blood flow rate by increasing them as needed based on an increase in patients metabolism. Failure to adjust the gas supply and the blood flow rate appropriately may cause insufficient O₂ supply needed or the amount of the patient's gaseous metabolism.
- Use of a pre-bypass filter is recommended to trap any particulate matter in the circuitry and priming solution.
- The Sampling Line is detachable. If the line is detached and not in use, please make sure to cap the end.
- Use of a recirculation line is recommended in the arterial line after the FX oxygenator to facilitate gaseous emboli removal during priming and emergency situations.
- It is recommended to use an air bubble detector between the reservoir outlet and FX oxygenator inlet. In the event gaseous emboli are detected, determine the source and take corrective action to eliminate further introduction of gaseous emboli into the perfusion circuit.
- Ensure the oxygenator is never subjected to a negative pressure. Negative pressure in the oxygenator can cause gaseous emboli to enter the blood phase from the gas phase.

Use CAPIOX FX25 Oxygenator with Hardshell Reservoir taking note of the additional precautions below.

- Use the safety devices, including the level sensor and the air detector, in order to control the volume of the blood in the reservoir and to detect gaseous emboli in the arterial lines.
- When rotating the venous blood inlet port on the header of Hardshell Reservoir, make sure that venous line and thermistor probe cable do not come into contact with any solution or administration line connected to the Hardshell Reservoir to prevent it from kinking.
- When separating Hardshell Reservoir from oxygenator module, hold the oxygenator firmly and remove connecting ring.
- The caps of unused ports should be left in place. This avoids contamination and prevents leakage of blood.
- Make sure the caps of unused luer ports are tightened firmly to prevent the leakage.
- The cardiotomy filter must be moistened in advance with the priming solution. If not moistened then this could compromise the effectiveness of the filter and the ability to use it up to the maximum blood flow rate.
- Where appropriate, introduce blood, plasma and/or blood derived products through quick prime port or any of the luer ports leading to cardiotomy filter after de-airing of oxygenator.
- Do not exceed the maximum operating volume of 4000 mL in the reservoir. This device is designed to be used with the volume of blood in the reservoir below 4000 mL.
- Avoid introduction of hemostatic agents into the Cardiotomy Reservoir/Oxygenator. Hemostatic agents are known to generate thrombi which may compromise performance of the Reservoir/Oxygenator.
- If blood entering the Cardiotomy Reservoir filter inlet appears to spill over versus entering the filter this may indicate that the filter is clogged. Discontinue using the Cardiotomy Reservoir filter. Replace the reservoir.
- Regularly check all tubing connections and ports for looseness or leakage.
- Scientific literature contains reports of altered dose response to some medications, such as nitroglycerin, fentanyl, etc., possibly due to varying degrees of absorption by synthetic materials such as those found in extracorporeal circuits. Avoid injection of medications into the Cardiotomy Reservoir filter.
- For a low dose of drug administration, do not inject from luer ports leading to the cardiotomy filter. If drug is injected inside the cardiotomy filter, it could stagnate in the filter.
- Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

Rx ONLY CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

METHOD OF OPERATION

Read WARNINGS (page 6) and PRECAUTIONS (page 7) before use. The following is a basic description of the use of the CAPIOX FX25 Oxygenator and Hardshell Reservoir.

a. SET-UP

1. Remove CAPIOX FX25 from package and examine it for defects.

CAUTION Do not use if the package or device is damaged (e.g. cracked) or any of the port caps are off.

2. Place CAPIOX FX25 in its holder.
 - a. CAPIOX FX25 Oxygenator with Hardshell Reservoir:
Set the holder at the appropriate height to prevent the arterial line and the blood cardioplegia line from kinking.
Making sure that the connecting ring is in place, set oxygenator in its holder (Code No.: XX*CXH18R, 801804 (Long Arm version) or 801139) with lever up, then push the lever down (See Fig. 1). If the oxygenator is not set into the holder correctly, the amount of fluid in the reservoir cannot be correctly assessed.
When rotating the oxygenator module, be sure to hold the support arm. If the oxygenator module is rotated while holding the oxygenator body, it could be damaged (See Fig. 2).

Fig. 1



Fig. 2

NOTE: When using CAPIOX FX25 Oxygenator module and Hardshell Reservoir separately,

- Hold the module firmly and separate it from Hardshell Reservoir by removing connecting ring.
- Set reservoir in the holder (Code No.:XX=CXH18R, 801804 or 801139) for CAPIOX FX25 Oxygenator with Hardshell Reservoir and pull the lever down. (see Fig. 1)
- Set the oxygenator module in its holder (Code No.:XX=CXH25F) and close the "lock" cover. (see Fig. 3-②)
- CAPIOX FX25 Oxygenator module without Hardshell Reservoir: Set the holder at the appropriate height to prevent the arterial line and the blood cardioplegia line from kinking. Set oxygenator module in its holder (Code No.:XX=CXH15) (See Fig.3-⑤)

CAUTION If the product is dropped during set-up, do not use it. Replace with another device.

- Set the sampling system into the holder.
- Connect water lines (1/2" (12.7 mm) tubing or couplers) to the water ports of CAPIOX FX25. (See Fig. 4-(1))

CAUTION Use the upper port for water supply and the lower port for water drainage. Otherwise heat exchange will not be performed sufficiently.

- Start water circulation through heat exchanger and circulate for at least 5 minutes. Check for leaks.

WARNING Do not use an oxygenator that leaks.

- Attach 1/2" (12.7 mm) venous line to the venous blood inlet port of venous reservoir. (See Fig. 4-(2))
- Attach 3/8" (9.5 mm) arterial line to the blood outlet port of CAPIOX FX25 Oxygenator. (See Fig. 4-(3))
- Remove the cap from the blood cardioplegia port and attach 1/4" (6.4 mm) blood cardioplegia line tubing. (See Fig. 4-(4)).

WARNING If the port is not to be used, push cap in and twist it to seal tightly or attach a short 1/4" (6.4 mm) tubing and clamp for security.

- Connect 3/8" (9.5 mm) pump line to the blood outlet port of Hardshell Reservoir. Connect other end to 3/8" (9.5mm) blood inlet port of oxygenator. (See Fig. 4-(5))
- The blood cardioplegia line must be connected to the luer port lateral to the blood outlet port of the oxygenator if oxygenated blood is needed for blood cardioplegia (see Fig. 4-(4)).
- Band all connections in the circuit.
- Connect 1/4" (6.4 mm) gas line to gas inlet port. (See Fig. 4-(6))

WARNING DO NOT OBSTRUCT GAS OUTLET PORT. (See Fig. 4-(7))

- Thermistor probes can be used with Y.S.I. 400 series* temperature monitor by use of the following parts.
Blue Cable (venous line): Code No.: CX+BP021
Red Cable (arterial line): Code No.: CX+BP022
* Y.S.I. 400 (trade name: Measurement specialties Inc.)
- Connect suction lines and vent line to suction ports on Hardshell Reservoir. Suction ports of CAPIOX FX25 Hardshell Reservoir have blue caps. Remove the blue caps before connection. (See Fig. 5)
- For use of CAPIOX FX25 Oxygenator module without Hardshell Reservoir, connect venous side male luer lock of sampling system to venous line. (See Fig. 6)
- When the luer port located on venous blood inlet port of CAPIOX FX25 Hardshell Reservoir is to be used, attach a three-way stopcock before use.
- Check figures and labels to make sure tubing is connected to the correct port. Special care should be taken to connect the correct tubing to the 1/4" (6.4 mm) gas inlet port (Fig. 4-(6)), and 1/4" (6.4 mm) blood cardioplegia port (Fig. 4-(4)).

CAUTION Ensure that all connected parts including the luer caps, the lock adapters and the port caps are securely affixed. Loose connections may cause contamination or a blood leak.

Fig. 3


Fig. 4


Fig. 5



Fig. 6

b. PRIMING PROCEDURE

CAUTION Perform priming using crystalloid solution which contains no blood, plasma and/or blood derived products. Using blood derived products will increase the amount of time necessary to prime. Therefore, if blood derived products are used during priming, start bypass after performing air purge thoroughly.

NOTE: CAPIOX FX25 can be primed without CO₂ flush when crystalloid priming solution is used. Use of CO₂ flush permits debubbling in a shorter time.

1. If recirculation line is used to prime, clamp arterial and venous lines distal to recirculation line, and make sure that recirculation line is not clamped.
2. Introduce crystalloid priming solution through quick prime port or any of the luer ports leading to cardiotomy filter.
3. Make sure the recirculation circuit and the purge line are not clamped, then start pump at a low speed. After checking for leakage or any other problem, gradually increase to full flow. Do not exceed 7 L/min flow rate. Vigorously recirculate the priming fluid through the entire circuit until all air bubbles are eliminated. Check oxygenator and tubing for leakage or any other problem. After all air bubbles are eliminated, circulate at full flow for 10 min to check oxygenator and tubing for leakage or any other problem.

WARNINGS

- Do not use an oxygenator and reservoir that leaks. Replace it with another CAPIOX FX25 oxygenator and reservoir.
- Do not use a tubing with internal diameter less than 3/16" (4.8 mm) as a recirculation line. Also do not use a sampling line nor a purge line for recirculation. If used, the oxygenator module can be damaged as a result of excessive positive pressure being generated inside itself.


CAUTIONS

- Do not supply gas during priming.
- Recirculate the priming solution at a rate of 4 L/min or higher to facilitate air removal. Failure to remove air from the oxygenator may result in serious injury to the patient.
- Maintain a minimum operating level of 200 mL in the reservoir at all times.
- Returning the priming solution back to the cardiotomy filter when there is an insufficient level of solution contained within the reservoir may generate gaseous emboli. Maintain an adequate level of solution in the reservoir.
- 4. When appropriate after the debubbling, introduce blood or blood-derived products through quick prime port or any of the luer ports leading to cardiotomy filter.
- 5. Set stopcocks as shown in Fig. 7, and close sampling line with arterial side stopcock to prevent arterial to venous shunting during extracorporeal circulation. After closing the purge line, reduce blood flow rate gradually to zero; then close recirculation line.
- 6. After priming, if air bubbles continue to appear, identify the cause and make necessary corrections. Remove the air while opening the purge line.

WARNING During recirculation, do not use pulsatile flow or stop blood pump suddenly. Otherwise gaseous emboli may enter the blood phase from the gas phase due to inertia force.

CAUTION Close purge line before stopping bypass.

Fig. 7

c. INITIATION OF BYPASS

Check the following before initiating bypass.

CAUTION Ensure that the de-airing process is complete prior to initiating bypass. Repeat b. "priming procedure" to dispel air.

Initiate extracorporeal circulation using normal procedure, taking note of the following warnings.

WARNINGS

- Start gas supply only after blood circulation is initiated.
- Before starting gas supply, confirm again that gas outlet port is not obstructed. Such obstruction may lead to pressure build-up in the gas phase, allowing gaseous emboli to enter into the blood phase.
- Before initiating extracorporeal circulation, be sure to confirm that recirculation line and purge line are closed and sampling line is also closed with arterial side stopcock. Otherwise, opening arterial line will cause the blood to flow back into reservoir through sampling line because of the patient's blood pressure and the head height.
- Start gas supply with V/Q=1 and FiO₂=100%, then make adjustments based on blood gas measurements.

d. DURING PERFUSION

1. In order to collect proper blood samples, withdraw at least 10 mL of blood, then collect blood through sampling line. In the case of arterial blood sampling, blood can be collected after opening the stopcock for arterial-venous shunting through sampling line.

WARNING Collect blood only while pump is running, or blood-side pressure will decrease and air bubbles may result.

NOTE: For the use of the sampling system separated from Hardshell Reservoir, sampling manifold holder (Code No: XX=XH051) is available.

2. Measure blood gases and make necessary adjustments as follows.
 - a. Control PaO₂ by changing concentration of oxygen in ventilating gas using gas blender.
 - To decrease PaO₂, decrease FiO₂.
 - To increase PaO₂, increase FiO₂.
 - b. Control PaCO₂ by changing the total gas flow.
 - To decrease PaCO₂, increase total gas flow.
 - To increase PaCO₂, decrease total gas flow.

WARNING A phenomenon called wet lung may occur when water condensation occurs inside fibers of microporous membrane oxygenators with blood flowing exterior to the fibers. This may occur when oxygenators are used for a longer period of time. If water condensation and/or a decrease in PaO₂ and/or an increase in PaCO₂ is noted during extended oxygenator use, briefly increasing the gas flow rate may improve the performance. Increase gas flow rate, to 20 L/min for 10 seconds. DO NOT repeat this flushing technique, even if oxygenator performance is not improved.

CAUTIONS

- A minimum of 0.5 L/min oxygen gas flow is needed when blood is circulated. Less than 0.5 L/min oxygen gas flow may result in inadequate gas exchange.
- Prior to resuming bypass, set FiO₂ at 100% to ensure adequate oxygenation. Increased PCO₂ and decreased PO₂ in the patient's blood at the beginning of recirculation may not be able to be recovered without appropriate gas supply.

3. Adjust blood flow rate, making sure that the blood level in the reservoir is adequate, to prevent gaseous emboli passing to patient. (When using CAPIOX FX25 with Hardshell Reservoir, refer to Fig. 11, Page 10 "SAFETY MARGIN".)
4. For purging, open the purge line.

CAUTION Close purge line at termination of bypass.

5. Do not give any shocks to the product during perfusion.

e. TERMINATION OF BYPASS

Terminate extracorporeal circulation using normal procedure, taking note of the following warnings:

WARNINGS

- Before reducing blood flow rate, confirm that purge line is closed and the sampling line is also closed with arterial side stopcock.
- On stopping pump, stop gas flow at once.
- If circulation is to be resumed, recirculation at a low flow rate is recommended. During recirculation, check blood gas pressure. An excessive gas flow may cause low PaCO₂, alkalosis, or blood damage.

f. OXYGENATOR REPLACEMENT

Always have a spare CAPIOX FX25 oxygenator for replacement if needed.

1. Set up and prime spare oxygenator module as described in SET-UP section. (page 7)
2. When necessary, maintain a decreased patient temperature as prescribed by physician.
3. Double-clamp blood inlet and outlet tubing connected to new oxygenator, and cut the tubing between the two clamps. (See Fig. 8)
4. Stop circulation, double-clamp venous and arterial line connected to old oxygenator, and cut the lines between the two clamps. (See Fig. 8)
5. Replace old oxygenator with new one by connecting venous and arterial line to the connectors at blood inlet and outlet tubing of new oxygenator. (See Fig. 9)

CAUTIONS

- After replacement, open recirculation line for debubbling.
- Band and secure all connections in the circuit.
- 6. Start circulation at a low flow rate.
- 7. Connect gas lines removed from old oxygenator module to the new module and start gas supply.
- 8. Connect water lines removed from old oxygenator module to the new module and start water supply. Check for leaks.
- 9. Connect the purge line of the new oxygenator to the luer port leading to the cardiotomy filter.
- 10. Connect sampling line of new oxygenator to venous side connecting site of old reservoir. (See Fig. 10)

Fig. 8


Fig. 9




Fig. 10

SAFETY MARGIN





Fig. 11 Safety Margin





* Minimum operating blood level in CAPIOX FX25 Hardshell Reservoir.

The above graph indicates the time in which the blood volume in the reservoir decreases to minimum operating blood level of 200 mL at each flow rate. If the blood volume in reservoir is below 200 mL, bubbles may pass to the oxygenator module.

PERFORMANCE DATA

Dedicated holders (sold separately)

Vacuum Assisted Venous Drainage (VAVD)

INTENDED USE

The FX25 Hardshell Reservoir is also intended for use in vacuum assisted venous drainage procedures. The use of a controlled vacuum on the hardshell improves venous drainage during minimally invasive surgical techniques or regular bypass surgery.

WARNINGS

- For drug administration in the hardshell reservoir, be aware that negative pressure in this line could draw more drug than intended. To prevent non intended drug amounts to be drawn, it is recommended to include only the intended amounts in the syringe.
- Always allow for a vented reservoir when vacuum is not being applied. Vent by declamping the vacuum release line.
- Avoid obstruction or occlusion of vacuum line to prevent retrograde air to the patient.
- Do not exceed -20 kPa (-150 mmHg) when using vacuum source for assisted venous drainage technique to prevent hemolysis.
- Open any port at the top of the reservoir for releasing negative pressure before stopping the pump. Otherwise the blood may flow sharply from the patient.
- When stopping the pump or when the blood flow rate is low during VAVD, all A-V shunt lines (ex. sampling line, purge line, etc) must be closed to prevent drawing air into the blood side of the oxygenator from the fiber and reversing the blood flow to the reservoir from the arterial side of the patient.
- The roller pump must be occluded properly because it is more likely to draw air into the blood side of the oxygenator from the fiber when VAVD is performed.
- Proceed slowly when returning to atmospheric pressure (elimination of vacuum). A sudden change in pressure may create turbulence of the blood inside the reservoir.
- Do not open the auxiliary port on the venous reservoir during VAVD to prevent bubble formation inside the venous filter.
- Do not occlude the positive pressure relief valve when vacuum assisted venous drainage is used. Over pressurization of the venous reservoir could obstruct venous drainage, decreasing the blood level in the venous reservoir and allowing air to enter both the venous and arterial side of the circuit. This could result in gaseous emboli entering the blood phase.

When VAVD is performed in conjunction with a centrifugal pump take caution as below.

- The line between the oxygenator and the centrifugal pump must be clamped before stopping the pump. Not clamping off the arterial line may draw air into the blood side of the oxygenator from the fiber. An arterial line one-way valve is recommended between the oxygenator and the centrifugal pump.
- Caution must be taken because the relationship between blood flow rate and pump speed changes depending on the degree of negative pressure applied.

PRECAUTIONS

- Do not connect the negative pressurized circuit to any other ports than the reservoir vent port. The negative pressurized circuit may suck blood back in the line.
- Do not leave the ports that are not in use open.
- A moisture trap is required due to the volume of condensation generated.
- Use a sterile negative pressurized circuit and do not reuse it.
- A positive pressure-warning device should be equipped on the reservoir.
- A controlled vacuum regulator is required.
- A positive pressure relief valve is required.
- The use of a negative pressure manometer gauge on the hardshell reservoir and a negative pressure relief valve [opening range at -20 kPa (-150 mmHg)] are recommended.
- Use of vacuum may cause liquid level to appear higher than actual due to reservoir compliance. A liquid level of 230 mL on volume scale assures actual minimum volume of 200 mL at a negative pressure of -20 kPa (-150 mmHg).

INSTRUCTIONS FOR USE

NOTE : 1-10 refers to the Structure diagram with numbered components and A-I refers to the Required Equipment for use with vacuum assisted venous drainage with letter indications.

1. Prepare CPB circuit set-up as per standard procedure.
2. Verify that blue caps from suction inlet ports 6 and auxiliary port 2 of reservoir are pushed all the way to insure adequate seal.
3. Secure in place all yellow luer caps; all caps are non-vented 4, 8.
4. Install controlled vacuum regulator D on wall suction source.
5. Remove positive pressure relief valve and attach a 3-way stopcock to the non-filtered luer lock port, and connect the sterile manometer gauge line E.
6. Connect the negative manometer gauge F to the sterile manometer line E.
7. Also, to the above mentioned (in procedure 5) 3-way stopcock, connect the positive pressure relief valve H to the other end of the 3-way stopcock.
8. Attach the gas filter I to the tubing connected to the moisture trap B and the "Y" connector.
9. Attach sterile moisture trap B tube to the vent port 5 on the hardshell reservoir.

<Initiating Bypass>

1. Begin with regular venous gravity drainage; at this time, the vacuum release line G is unclamped.
2. To initiate vacuum assist venous drainage: set the vacuum regulator D to -5.3kPa (-40 mmHg), then clamp the vacuum release line G.
3. Monitor the negative pressure inside the hardshell reservoir with the negative manometer gauge F.
4. Adjust the negative pressure to optimize venous return. Set vacuum regulator D between -5.3 kPa and -8.0 kPa (-40 and -60 mmHg).

NOTE:

- The "DO NOT OBSTRUCT" caution label near the vent port is not applicable for vacuum assisted venous drainage.
- Use an occlusive roller head for the suction and LV vent lines.

<Coming off bypass>


Unclamp vacuum release line G, venous return will quickly fall. Come off bypass as per standard procedure.

REQUIRED EQUIPMENT


- A FX25R venous reservoir
- B Sterile moisture trap
- C Clamp
- D Controlled vacuum regulator [adjustable between 0 to -20 kPa (-150 mmHg)]

- E Sterile manometer gauge line
- F Negative manometer gauge
- G Sterile vacuum release line
- H Sterile positive pressure relief valve [opening range from 0 to 1.3 kPa (10 mmHg)]
- I Sterile gas filter

SET-UP CONFIGURATION

STRUCTURE

Post – Operative Chest Drainage

CONTRAINDICATIONS

Post-operative chest drainage with subsequent autotransfusion procedures are contraindicated in the following events:

- Gross perforations to the chest wall or an air leak in the lungs.
- Infection or malignancy occurring systemically or in the pericardium, mediastinum, or lungs.
- Suspected or evident gross contamination with foreign material, lymphatic failure, or perforated intestine.

- Presence of the following in the aspiration site: topical hemostatic agents, bactericidal wound irrigants or antibiotics not intended for parenteral administration.
- Open chest and vacuum applied.
- Administration of protamine prior to reservoir being removed from bypass circuit.
- Patient returned to surgery for any reason.
- Use of vented chest tubes that do not have vent flow regulation such as stopcocks.

WARNINGS

- A qualified person should assess the quality and suitability of returning any blood that has been collected before reinfusion begins. The safe reinfusion of collected fluids is the sole responsibility of the attending physician.
- Complications such as the following have been associated with chest drainage and subsequent reinfusion: blood trauma, blood coagulation, coagulopathies, and particulate or air embolism.
- Reinfusion of collected blood/fluid should be done on an hourly basis unless fewer than 50 mL are collected hourly.
- A minimum reservoir level of 20 mL of fluid should always be maintained to prevent air embolus passing to patient.
- Blood that has been in the reservoir for longer than 4 hours should not be transfused.¹
- It is not recommended that autotransfusion continue longer than 18 hours after surgery.²
- Attach a shunt line (the bridge) to connect the filtered and unfiltered sections in the reservoir as indicated in instructions for use.

- When using reservoir, monitor the attached shunt line for blood in the line. If blood is detected in the line this is an indication that the reservoir filter has clotted. Replace reservoir promptly.
- Filter occlusion during high volume chest drainage may cause blood/fluid to pass through the external shunt line bypassing filtration. All fluids that pass through the shunt line must be filtered prior to reinfusion.
- Accepted medical and nursing care routines must be followed during chest drainage.
- When vacuum is used during chest drainage, do not exceed -20 kPa (-150 mmHg) (-195 cmH₂O).

¹ Standards for Blood Banks and Transfusion Services, 16th ed. American Association of Blood Banks

² Page. R, et al, Hard-Shell Cardiomy Reservoir for Reinfusion of Shed Mediastinal Blood. Ann Thorac Surg 1989; 48:514-7.

Hardshell Reservoir Operating Specifications for Post-Operative Chest Drainage

ITEM	SPECIFICATION
Cardiotomy Maximum Blood Flow Range	5 L/min
Blood Storage Capacity	4,000 mL
Minimum level in reservoir	20 mL

List of Required Materials for Use of Reservoir for Post-Operative Chest Drainage

*REFERENCE	QUANTITY	DESCRIPTION	TUBING LENGTH
[a]	6 each	1/4" (6.4 mm) I.D. dead-end tubing with cap	2" (50.8 mm) each
[b]	1 each	3/8" (9.5 mm) I.D. tubing with caps	12" (304.8 mm)
[c]	1 each	1/2" (12.7 mm) X 1/2" (12.7 mm) connector with 1/2" (12.7 mm) I.D. dead-end tubing with cap	6" (152.4 mm)
[d]	2 each	3/8" (9.5 mm) X 3/8" (9.5 mm) Y connector with 3/8" (9.5 mm) I.D. 40 Durometer tubing with caps	96" (2438.4 mm)
[e]	1 each	3/8" (9.5 mm) I.D. dead-end tubing with cap	4" (101.6 mm)
[f]	1 each	3/8" (9.5 mm) perfusion adapter attached to 4" (101.6 mm) of 1/8" (3.175 mm) I.D. tubing that includes a clamp with a female luer end, connected to male luer with 48" (1219.2 mm) of 1/8" (3.175 mm) I.D. tubing with another female luer end connected to a male luer with 2.5" (12.7 mm) of 3/16" (4.7625 mm) I.D. tubing; the end of the tubing is covered by a cap.	N/A
[g.1.]	1 each	Pressure relief valve: opening vacuum range: 3.1 - 6.2 kPa (30 - 60 cmH ₂ O) (23.1 - 46.2 mmHg) with bacterial filter	N/A
[g.2.]	1 each	Pressure relief valve: opening pressure of approximately 0.8 kPa (7.9 cmH ₂ O) (6 mmHg)	N/A
[h]	1 each	1/4" (6.4 mm) I.D. tubing	As needed
[i]	1 each	Water manometer	N/A
[j]	1 each	1/4" (6.4 mm) I.D. tubing	18" (457.2 mm)
[k]	1 each	Male non-vented luer lock cap	N/A
[l]	1 each	Three way stopcock capable of opening three ways simultaneously	N/A

INSTRUCTIONS FOR USE

NOTE: When referring to the instructions for Post Operative Chest Drainage, the bracketed letters refer to the List of Required Materials on page 13 and the parenthesized numbers refer to the illustration on page 12.

- Carefully remove tubing and/or blue caps from suction inlet ports (6) of reservoir. Plug five of the suction inlet ports with 1/4" (6.4 mm) I.D. dead-end tubes [a].
- Attach 3/8" (9.5 mm) I.D. tubing [b] from the auxiliary port (2) to one of the suction inlet ports (6); all suction inlet ports are 1/4" (6.4 mm) - 3/8" (9.5 mm). This will create a bridge between the filtered and unfiltered sections of the reservoir.
- Clamp the 1/2" (12.7 mm) circuit tubing attached to the 1/2" (12.7 mm) venous blood inlet port (10) of the reservoir, cut it and connect a 1/2" (12.7 mm) connector with 1/2" (12.7 mm) I.D. dead-end tube [c].
- Connect the chest drainage tubes to the 3/8" (9.5 mm) X 3/8" (9.5 mm) Y connector with 96" (2438.4 mm) of 3/8" (9.5 mm) 40 Durometer tubing [d]; if only one chest tube is used, plug one of the 3/8" (9.5 mm) ports of the Y connector with a 3/8" (9.5 mm) dead-end tube [e]. Attach the other end of the 3/8" (9.5 mm) X 3/8" (9.5 mm) Y connector with 96" (2438.4 mm) of tubing to the vertical port to CR filter (7). In case 3 or more chest tubes are used, attach an additional 3/8" (9.5 mm) X 3/8" (9.5 mm) Y connector with 96" (2438.4 mm) of 3/8" (9.5 mm) 40 Durometer tubing [d] to one of the suction inlet ports (6) that was capped off by a dead-end tube.
- Secure in place all yellow luer caps; all caps are non-vented (4, 8).
- Disconnect sampling system (3) from the back of the venous blood inlet port (10) and attach a non-vented luer lock cap [k].
- Attach a 1/4" (6.4 mm) dead-end tube [a] to the quick prime port (9) on the reservoir.
- Hold the reservoir parallel to the floor with the venous blood outlet port at the top to prevent blood spillage. Clamp the circuit tubing attached to the outlet port, cut it, then connect a 3/8" (9.5 mm) to 1/8" (3.175 mm) blood infusion adapter line [f]. Clamp the adapter line and leave the end cap in place until autotransfusion is set up. In case of using 1/4" (6.4 mm) adapter, take 1/4" (6.4 mm) adapter off first.
- Remove positive pressure relief valve from non-filtered luer lock port and attach three-way stopcock. Secure to the stopcock the negative pressure relief valve with attached bacterial filter [g.1] and positive pressure relief valve [g.2].
- Install the hardshell reservoir into the autotransfusion holder (product code: XX*XH032).

CAUTION Maintain the hardshell reservoir below the thorax to facilitate chest drainage.

- Using a 1/4" (6.4 mm) I.D. tubing [h], connect the vent port (5) on top of the hardshell reservoir to a regulated vacuum source with 1.5 to 2.1 kPa (15-20 cmH₂O) (11.6-15.4 mmHg).

NOTE: The "DO NOT OBSTRUCT" caution label near the vent port is not applicable for post-operative chest drainage and autotransfusion.

OR

- The use of a water seal/water manometer [i] is highly recommended for chest drainage. If using a water manometer, connect one end of a 1/4" (6.4 mm) I.D. tubing [j] to the vent port (5) and the other end to the water manometer.

WARNING To minimize blood damage, vacuum should be regulated so as not to exceed 13.3-20 kPa (100-150 mmHg) (130-195 cmH₂O).

12. Record the time drainage into the reservoir begins.

WARNING It is recommended that autotransfusion continue no longer than 18 hours after surgery.

13. Continually monitor drainage and elapsed time. Determine hourly drainage rate.

WARNING Reinfusion of collected blood/fluid should be done on an hourly basis unless fewer than 50 mL are collected hourly.

14. Upon collection of 50 mL of blood, remove the end cap of the blood infusion adapter line and attach the blood infusion adapter to an I.V. Administration set/pump apparatus.

NOTE: It is recommended that the infusion pump be used with an air detection sensor for all reinfusion procedures.

15. Uncamp the blood infusion line and proceed with autotransfusion.

WARNINGS

- To reduce risk of air embolism, evacuate all air completely from the infusion line prior to autotransfusion.
- To minimize blood damage, do not exceed 13.3-20 kPa (100-150 mmHg) (130-195 cmH₂O) pressure when reinfusing blood.
- The physician is responsible for determining the optimal time, rate and volume to be autotransfused to the patient.
- The clinician should be aware of potential complications such as air emboli, associated with the use of the I.V. Administration set and the infusion pump. Use the I.V. Administration set and the infusion pump in accordance with the instructions provided by the manufacturers.

SPECIFICATIONS

Product Description

Reorder Number	Description
CX*FX25RW	Oxygenator with integrated arterial filter with Hardshell Reservoir
CX*FX25RE	Oxygenator with integrated arterial filter with Hardshell Reservoir
CX*FX25W	Oxygenator with integrated arterial filter
CX*FX25E	Oxygenator with integrated arterial filter

(Refer to "DIRECTION OF BLOOD OUTLET PORT", page 5)

CAPiox FX25 Oxygenator Module

COMPONENT	SPECIFICATIONS	
Housing	Material	Polycarbonate
Fibers	Material	Microporous Polypropylene
	Surface area	Approx. 2.5 m ²
Arterial filter	Material	Polyester screen type
		Pore size 32 µm
Heat Exchanger	Material	Stainless steel
	Surface area	Approx. 0.2 m ²
Blood Flow Range		Min. 0.5 L/min Max. 7.0 L/min
Priming Volume (Static)		260 mL
Blood inlet port (from pump)		3/8" (9.5 mm)
Blood outlet port		3/8" (9.5 mm)
Cardioplegia port		1/4" (6.4 mm)
Gas inlet port		1/4" (6.4 mm)
Gas outlet port		1/4" (6.4 mm)
Water ports		1/2" (12.7 mm) Hansen quick connect fittings
Maximum Pressure	Blood inlet	133 kPa (1,000 mmHg)
	Water inlet	196 kPa (2 kgf/cm ²)

Hardshell Reservoir Section

COMPONENT	SPECIFICATIONS	
Housing	Material	Polycarbonate
Blood Flow Range		Venous Flow : Min. 0.5 L/min Max. 7.0 L/min Cardiotomy Inlet : Max. 5.0 L/min Combined Flow : Max. 7.0 L/min
Blood Storage Capacity		4,000 mL
Minimum Operating Volume		200 mL
Venous Filter	Material	Polyester screen type
		Pore size 47 µm
Cardiotomy Filter	Material	Polyester depth type
Defoaming Part	Material	Polyurethane foam
Venous blood inlet port		1/2" (12.7 mm) rotatable
Blood outlet port (to pump)		3/8" (9.5 mm)
Suction ports		six 1/4" (6.4 mm)
Vertical port to CR filter		3/8" (9.5 mm)
Quick prime port		1/4" (6.4 mm)
Vent port		1/4" (6.4 mm)
Auxiliary port		1/4"—3/8" (6.4 mm—9.5 mm)
Luer ports		-three filtered luer locks to cardiotomy filter -a non-filtered luer lock -two luer locks on venous inlet
Maximum Sustainable Negative Pressure in Reservoir		-20 kPa (-150 mmHg)

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LIMITED WARRANTY

TERUMO Corporation warrants that the product has been manufactured and packaged with reasonable care, and will be free from defects under normal and proper use. TERUMO Corporation will not be liable for any direct, indirect, incidental, special, general or consequential loss, damage, or expense from the use of this product. TERUMO Corporation's sole obligation under this warranty shall be to repair or replace, at its option, any products that it believes were defective at the time of shipment if notice thereof is received before the expiration date indicated on the package. THIS WARRANTY IS THE SOLE AND ONLY WARRANTY OF TERUMO CORPORATION AND IS EXPRESSLY IN LIEU OF AND EXCLUDES ANY OTHER EXPRESSED OR IMPLIED WARRANTY OF ANY KIND, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. IN ADDITION, THE REMEDIES SET FORTH HEREIN ARE THE SOLE REMEDIES AVAILABLE TO ANY PERSON OR ENTITY. No agent, employee or representative of TERUMO Corporation has any authority to change, amend or alter the foregoing or assume or bind TERUMO Corporation to any additional liability or responsibility in connection with this product.

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©TERUMO CORPORATION 2013-07
Xcoating is a trademark of TERUMO CORPORATION.