


PLAJEX Syringe I

Do not reuse

【CAUTION】
Do not reuse or re-sterilize

【DESCRIPTION】

<Component>



<Materials which attach to inside the body>

Component	Material
Syringe barrel	Cyclo olefin polymer
Stopper	Chloro butyl rubber
Needle	Stainless steel, The mixture of silicone oil

【INTENDED USE】

PLAJEX Syringe I is a syringe barrel for use as a medical drug container designated as “ready-to-be filled”.
Terumo is manufacturing empty pre-fillable syringe components only which will be further processed by Pharmaceutical companies as per their filling requirements.

【INSTRUCTIONS FOR USE】

1. Pull on gloves if necessary.
2. Take a PLAJEX Syringe I out of the package.
3. Disinfect the injection site.
4. Remove the needle shield and the rigid needle shield by pulling straight without twisting.
5. Insert the needle into the injection site and deliver the drug.
6. Pull the needle out and cover the injection site with a small bandage if necessary.

<Precautions for use>

- (1) While pulling off the needle shield (and the rigid needle shield), do not let the needle point touch any surface to avoid the needle point being deformed or blunt.
- (2) Do not touch the needle directly with your hands to prevent from accidental needle-stick injury or infection.

【PRECAUTIONS】

<Important precautions>

- (1) Do not apply excessive force on the needle or the joint between the needle and the syringe barrel. [The excessive force may cause damage to the needle or the joint between the needle and the syringe barrel, which can lead to liquid leakage.]
- (2) When you re-cap the needle after use, use a protection device or re-cap by placing the shield(s) on a stand, not holding it with hands, and directing the needle straight into the shield to prevent from accidental needle-stick injury. [Inclined needle shield (and Rigid needle shield) may cause the needle to penetrate though the needle shield (or the rigid needle shield)]
- (3) Make sure to pull the plunger rod straight. [If you pull the plunger rod at a slant, it may decrease adhesion between the stopper and the barrel, which can lead to liquid leakage, air introduction, or detached stopper.]
- (4) Do not apply such force as strongly gripping the barrel or other part, or compressing or bending the syringe. [Such force may decrease adhesion between the stopper and the barrel, which can lead to liquid leakage, air introduction, or damages.]
- (5) Do not move the plunger rod repeatedly back and forth.
- (6) When storing PLAJEX in a refrigerator, handle with care. [In low temperature the shock resistance strength may be decreased, which can lead to breakage of the syringe.]
- (7) Do not reuse the same syringe and do not administer the residual liquid.

【STORAGE AND EXPIRATION】

<Storage>

Keep PLAJEX Syringe I dry and avoid direct sunlight, excessive heat and humidity.

<Expiration>

Expiration is printed on the carton box. (By self certification)

【CONTACT】

Terumo Corporation
Telephone +81-0120-12-8195 Terumo Call Center

